


Medical Tape

[Non-woven]



COMPOSITION:


Non-woven fabric,Acrylic adhesive,Backing paper

BACKING PAPER CUTTING TYPE:


Line slitting/S-shaped slitting

APPLICATION:

Used to provide adhesive force for wound dressings, bandages, infusion tubes, etc.



Different wound parts of the body



Width	Length	Box size	Packaging	Customization
5cm	10m	5.5x7x7cm	150roll/carton	Size Shape Material
10cm	10m	13.5x7x7cm	100roll/carton	
15cm	10m	18x7x7cm	100roll/carton	
20cm	10m	20.5x7x7cm	100roll/carton	

