

Custom PU Dressings

Your brand. Your format. Your peel design.



Rolls. Sheets. Custom shapes.

A specification-led approach to printed PU film dressings, from material and artwork review to custom sample preparation.

Taizhou Longer Intelligent Technology Co., Ltd.

Product images illustrate customization concepts, not fixed stock specifications.

Built around your dimensions.

Examples from customer projects. These are not a fixed size range.



Rolls

10 cm x 5 m

Most frequently supplied size. A 20 cm x 5 m version has also been sampled.



Sheets

20 x 30 cm | 20.8 x 24.8 cm

Customer-specific dimensions. Other sizes can be developed to the project brief.



Shaped films

Oval: 15 x 10 cm

Outline, central opening and peel layout are reviewed against the customer drawing.

Each layer has a role.

A typical construction, with options reviewed for the requested design.

01

Surface carrier

PET or PP; gloss or matte finish.

02

PU film

Generally kept to the established specification.

03

Acrylic adhesive

Coat weight is generally kept fixed.

04

Glassine release liner

80 g/m² is typical; 100 g/m² or other weights on request.

Schematic only - layer thicknesses are not to scale. Final material selection is confirmed during design review.

Designed on both sides.

Artwork and cut layouts are developed together for the agreed construction.



Front printing

Branding and artwork can be specified on the front surface.

Back printing

Back-side printing can also be requested. Feasibility depends on the selected construction and requires confirmation.

Front peel cut

A central S-shaped cut is the usual design. Alternative cut layouts can be reviewed.

Back peel cut

A straight edge kiss-cut is the usual design. Position and geometry can be customized.

From your brief to a sample.

We review your requirements and provide practical design recommendations.

01

Share the requirement

Send the format, dimensions, drawings or reference sample, artwork and expected quantity.

02

Review the design

Confirm PET or PP, finish, liner weight, print feasibility and the front and back peel layouts.

03

Prepare the sample

Custom sample preparation is available. Timing is confirmed after the design review.

04

Confirm for supply

Review the sample and agree the final specification, quantity and supply arrangements.

Flexible MOQ

Discussed according to the project requirements.

Tell us what you need.

A drawing, reference sample or clear brief is a useful starting point.

Your enquiry checklist

- Format and dimensions
- Surface carrier and finish
- Front and back printing requirements
- Front S-cut and back edge kiss-cut requirements
- Release liner weight and expected quantity

Supply basis

Generally supplied non-sterile. Customers typically arrange subsequent packaging and sterilization. Final requirements are agreed for each project.

Talk to Longer iTech

WhatsApp +86 133 8213 3835

Taizhou Longer Intelligent Technology Co., Ltd.

Building 7, North Side of Kechuang Pioneering Avenue,
Gaogang District, Taizhou City, Jiangsu Province, China