


Wound dressing

[Non-woven + PU]




COMPOSITION:

Non-woven fabric/PU film,Acrylic adhesive, Cotton pad,Release paper,Protection paper

APPLICATION:

- 1.Protection for minor scratches, small area first and second degree burns, and minor wounds.
- 2.Prevention and treatment of bedsores and skin cracking.
- 3.Used for protection of donor site wounds and surgical incisions.



6x7cm




9.5x6.5cm




10x10cm



10x10cm




10x15cm



10x20cm



10x10cm



10x35cm